

# HVO Pro Gen 10 Oxygen Concentrator



## Specifications

Flow Rate	Up to 10 LPM
Dimensions	14.6" W x 13.4" D x 24" H / 372 x 340 x 612mm
Electrical Requirements	120V/60Hz (± 10%)
Weight	53 lbs / 24kg
Oxygen Concentration	93% ± 3 (V/V)
Operating Temperature	50 to 90 F / 12.8 to 32.2 C
Outlet Pressure	5.8-11.6 psi / 40-80 Kpa
Standard Features	Low Noise Design HEPA filter Flame retardant function Suitable for 24x7 operation
Warranty	1 year
Storage/Transport Temperature	-30 to 160 °F
Operating Temperature	55 to 90 °F
Storage/Transport Humidity	Up to 95%, non-condensing
Visual Alarm Indicators	High system pressure Low system pressure Power failure