



O2 CONCENTRATOR

QUALITY

- ✓ Triple air filtration including carbon and bacterial filtration.
- ✓ Multi parameter display, including percentage of concentration.
- ✓ Easy filter access.
- ✓ Consistent supply of high concentration supplemental oxygen.

ECONOMICAL

- ✓ Has the ability to pay for itself in less than a year.
- ✓ Low total cost of ownership.
- ✓ Eliminates the requirement of delivered oxygen.

SAFE

- ✓ Power failure, no flow, high and low pressure alarms.

CONVENIENT

- ✓ Smooth quiet operation.
- ✓ Casters for portability.
- ✓ Programable digital run timer.
- ✓ Easily adaptable to any anesthesia system (O2 DISS MALE and Hose Barb fitting included)

SPECIFICATIONS

Weight	5 lpm: 55 lbs OR 10 lpm: 60 lbs
Height	25 in
Depth	12.5 in
Width	15 in
Flow Meter Capacity	0-5 lpm OR 0-10 lpm
Sound	42 - 50 dBa
Line Pressure	10 psi



High Volume Oxygen
 500 W South Street, Ste #4
 Lincoln, Nebraska 68522

Phone: (402) 476-0555

Email: sales@highvolumeoxygen.com